## **Product Data Sheet**

к

## PerCP/Cy5.5 anti-human CD26

| Catalog # / Size:      | 2113580 / 100 tests<br>2113575 / 25 tests   |   |
|------------------------|---|---|
| Clone:                 | BA5b  | 27 M  |
| Isotype:               | Mouse IgG2a, к  |   |
| Reactivity:            | Human   |   |
| Preparation:           | The antibody was purified by affinity<br>chromatography and conjugated with<br>PerCP/Cy5.5 under optimal conditions.<br>The solution is free of unconjugated<br>PerCP/Cy5.5 and unconjugated<br>antibody. | Relative Cell   |
| Formulation:           | Phosphate-buffered solution, pH 7.2,<br>containing 0.09% sodium azide and<br>0.2% (w/v) BSA (origin USA).   | Log Fluorescence Intensity<br>Human peripheral blood<br>lymphocytes were stained with |
| Workshop<br>Number:    | VI N-L078   | CD26 (clone BA5b) PerCP/Cy5.5<br>(filled histogram) or mouse IgG2a,                   |
| <b>Concentration</b> : | Lot-specific  | PerCP/Cy5.5 isotype control (open<br>histogram).                                      |

## **Applications:**

| Applications: | Flow Cytometry |
|---------------|----------------|
|---------------|----------------|

**Recommended** Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

Application
References:
Schlossman S, *et al.* Eds. 1997. Leucocyte Typing VI. Garland Press. London.
Schlossman S, *et al.* Eds. 1995. Leucocyte Typing V. Oxford University Press. New York.

**Description:** CD26 is a 110 kD type II membrane protein also known as ADA-binding protein and dipeptidyl peptidase IV (DPPIV). It is a member of the peptidase and ectoenzyme family. CD26 is expressed on the membrane of mature thymocytes, T lymphocytes (upregulated upon activation), B cells, NK cells, and macrophages. CD26 cleaves off N-terminal X-Pro and X-Ala dipeptides from polypeptides. It plays an integral role as a costimulatory molecule in T cell activation. CD26 may interact with extracellular matrix proteins such as fibronectin or collagen, CD45 and ADA.

| Antigen            | 1. Kameoka J, <i>et al.</i> 1993. <i>Science</i> 261:466.   |
|--------------------|---|
| <b>References:</b> | 2. Dang N, <i>et al.</i> 1990. <i>J. Exp. Med.</i> 172:649. |

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com