## **Product Data Sheet**

## Pacific Blue<sup>™</sup> anti-human CD45

Catalog # / Size:	2120145 / 100 tests 2120105 / 25 μg 2120110 / 100 μg	
Clone:	НІЗО	La C
		Env.
Isotype:	Mouse IgG1, к	10
<b>Reactivity:</b>	Human	Relative Cell Number
Preparation:	The antibody was purified by affinity chromatography, and conjugated with Pacific Blue™ under optimal conditions. The solution is free of unconjugated Pacific Blue™.	10 <sup>0</sup> 10 <sup>1</sup> Log Flu
Formulation:	test size: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin USA). microg size: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.	Human periph lymphocytes s Pacific Blue™
Workshop Number:	IV N816	
Concentration:	test size: lot-specific; microg sizes: 0.5 mg/ml	

10<sup>0</sup> 10<sup>1</sup> 10<sup>2</sup> 10<sup>3</sup> 10<sup>4</sup> Log Fluoresence Intensity

Human peripheral blood lymphocytes stained with HI30 Pacific Blue™

## **Applications:**

Flow Cytometry
Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. <b>For test size</b> , the suggested use of this reagent for immunofluorescent staining is 5 microL per $10^6$ cells in 100 microL volume. <b>For microg sizes</b> , the suggested use of this reagent for immunofluorescent staining is $\leq 0.5$ microg per $10^6$ cells in 100 microL volume. It is recommended that the reagent be titrated for optimal performance for each application.
* Pacific Blue™ has a maximum emission of 455 nm when it is excited at 405 nm. Prior to using Pacific Blue™ conjugate for flow cytometric analysis, please verify your flow cytometer's capability of exciting and detecting the fluorochrome.
Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen tissue sections and formalin-fixed paraffin-embedded tissue sections <sup>9</sup> , inhibition of CD45 functions4, immunofluorescence <sup>11</sup> , and Western blotting3.
It was found that the HI30 clone and the 2D1 clone can cross block each other's binding.
<ol> <li>Knapp W, <i>et al.</i> 1989. Leucocyte Typing IV. Oxford University Press. New York.</li> <li>Kishihara K, <i>et al.</i> 1993. <i>Cell</i> 74:143.</li> <li>Esser M, <i>et al.</i> 2001. <i>J. Virol.</i> 75:6173. (WB)</li> <li>Yamada T, <i>et al.</i> 2002. <i>J. Biol. Chem.</i> 277:28830.</li> <li>Nagano M, <i>et al.</i> 2007. <i>Blood</i> 110:151.</li> </ol>

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com

6. Jiang Q, <i>et al.</i> 2008. <i>Blood</i> 112:2858. <u>PubMed</u> 7. Morozov A, <i>et al.</i> 2010. <i>Clin Cancer Res.</i> 16:5630. <u>PubMed</u>
8. Yoshino N, <i>et al.</i> 2000. <i>Exp. Anim. (Tokyo)</i> 49:97. (FC)
9. Friedman T, <i>et al.</i> 1999. <i>J. Immunol.</i> 162:5256. (IHC)
10. Oeztuerk-Winder F, <i>et al.</i> 2012. <i>EMBO J.</i> 31:3431. (FC) <u>PubMed</u>
11. Rees LE, <i>et al.</i> 2003. <i>Clin. Exp. Immunol.</i> 134:497. (IF)
12. Lee J, <i>et al.</i> 2015. <i>J Exp Med.</i> 212:385. <u>PubMed</u>
13. Breton G, <i>et al.</i> 2015. <i>J Exp Med.</i> 212:401. <u>PubMed</u>
14. Marquardt N, <i>et al.</i> 2015. <i>J Immunol.</i> 6:2467. <u>PubMed</u>

**Description:** CD45 is a 180-240 kD single chain type I membrane glycoprotein also known as leukocyte common antigen (LCA) and T200. It is a tyrosine phosphatase expressed on the plasma membrane of all hematopoietic cells, except erythrocytes and platelets. CD45 is a signaling molecule that regulates a variety of cellular processes including cell growth, differentiation, cell cycle, and oncogenic transformation. CD45 plays a critical role in T and B cell antigen receptor-mediated activation by dephosphorylating substrates including p56Lck, p59Fyn, and other Src family kinases. CD45 non-covalently associates with lymphocyte phosphatase-associated phosphoprotein (LPAP) on T and B lymphocytes. CD45 has been reported to bind galectin-1 and to be associated with several other cell surface antigens including CD1, CD2, CD3, and CD4.

Antigen	1. Thomas M. 1989. <i>Annu. Rev. Immunol.</i> 7:339.
<b>References:</b>	2. Trowbridge I, et al. 1994. Annu. Rev. Immunol.12:85.