Product Data Sheet

Brilliant Violet 605[™] anti-human CD62P (P-Selectin)

Catalog # / Size:	2124600 / 100 tests 2124595 / 25 tests	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
Clone:	AK4	Interpret to the second state of the second
Isotype:	Mouse IgG1, κ	
Reactivity:	Human	
Preparation:	The antibody was purified by affinity chromatography and conjugated with Brilliant Violet 605 [™] under optimal conditions. The solution is free of unconjugated Brilliant Violet 605 [™] and unconjugated antibody.	
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA).	
Workshop Number:	VI P-44	
Concentration:	Lot-specific	

Applications:

Applications:	Flow Cytometry	
Recommended Usage:	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤ 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.	
	Brilliant Violet 605 [™] excites at 405 nm and emits at 603 nm. The bandpass filter 610/20 nm is recommended for detection, although filter optimization may be required depending on other fluorophores used. Be sure to verify that your cytometer configuration and software setup are appropriate for detecting this channel. Refer to your instrument manual or manufacturer for support. Brilliant Violet 605 [™] is a trademark of Sirigen Group Ltd.	
	This product is subject to proprietary rights of Sirigen Inc. and is made and sold under license from Sirigen Inc. The purchase of this product conveys to the buyer a non-transferable right to use the purchased product for research purposes only. This product may not be resold or incorporated in any manner into another product for resale. Any use for therapeutics or diagnostics is strictly prohibited. This product is covered by U.S. Patent(s), pending patent applications and foreign equivalents.	
Application Notes:	Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen tissue sections4 and <i>in</i> <i>vitro</i> blocking of adhesion of platelets1. The LEAF ^{m} purified antibody (Endotoxin <0.1 EU/µg, Azide-Free, 0.2 µm filtered) is recommended for functional assays (Cat. No. 304911).	
Application References:	1. Skinner M, <i>et al.</i> 1991. <i>J. Biol. Chem.</i> 266:5371. (Block) 2. Kishimoto T, <i>et al.</i> Eds. 1997. Leucocyte Typing VI. Garland Publishing Inc. London.	

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com 3. Yen YT, *et al.* 2006. *J. Virol.* 80:2684. 4. Sato Y, *et al.* 2005. *Blood* 106:428. (IHC)

Description: CD62P is a 140 kD type I transmembrane glycoprotein also known as P-selectin, platelet activation-dependent granule membrane protein (PADGEM), and GMP-140. It is expressed on activated platelets, megakaryocytes, and endothelial cells. CD62P is primarily stored in secretory α-granules in platelets and Weibel-Palade bodies in endothelial cells, and is rapidly relocated to the plasma membrane upon activation. The ligands for CD62P are CD162 and CD24. A primary function of CD62P is cell adhesion during neutrophil rolling, and platelet-neutrophil and platelet-monocyte interactions.

 Antigen
 1. McEver R, et al. 1995. J. Biol. Chem. 270:11025.

 References:
 2. Varki A. 1994. P. Natl. Acad. Sci. USA 91:7390.