Alexa Fluor[®] 647 anti-human CD162

Catalog # / Size:	2244050 / 100 tests 2244045 / 25 tests	
Clone:	KPL-1	selative Cell Number
lsotype:	Mouse IgG1, к	
Immunogen:	PSGL-1 transfected murine 300.19 pre B-cell line	
Reactivity:	Human, Non-human primate, Other	
Preparation:	The antibody was purified by affinity chromatography and conjugated with Alexa Fluor® 647 under optimal conditions. The solution is free of unconjugated Alexa Fluor® 647.	ř.
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin USA).	Human stained Alexa Fl histogra isotype Alexa Fl
Concentration:	Lot-specific	

Human peripheral blood was stained with CD162 (clone KPL-1) Alexa Fluor® 647 (filled histogram) or mouse IgG1, ĸ isotype control (clone MOPC-21) Alexa Fluor® 647 (open histogram).

0 10³ 10⁴ Log Fluorescence Intensity

105

Applications:

Applications:	Flow Cytometry
Recommended Usage:	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 µl per million cells or 5 µl per 100 µl of whole blood. It is recommended that the reagent be titrated for optimal performance for each application. * Alexa Fluor® 647 has a maximum emission of 668 nm when it is excited at 633 nm / 635 nm.
Application Notes:	Clone KPL-1 is reported to recognize the tyrosine sulfation consensus motif of PSGL-1 ¹ . Additional reported applications (for the relevant formats) include: Western Blot ¹ , immunoprecipitation ² , immunohistochemical staining of acetone-fixed frozen tissue sections and formalin-fixed paraffin embedded tissue sections ¹ , blocks the recognition of PSGL-1 with P- and L- selectin ¹ .
Application References:	1. Snapp KR, <i>et al.</i> 1998. <i>Blood</i> 91:154.
Description:	CD162, also known as p-selectin glycoprotein ligand-1 (PSGL-1), is a 120 - 220 kD, mucin-like type I transmembrane glycoprotein. CD162 binds to CD62P (P-Selectin), CD62E (E-Selectin) and CD62L (L-Selectin). The interactions between P-selectin and P-selectin glycoprotein ligand-1 (PSGL- 1) mediate the earliest "rolling" of leukocytes on the lumenal surface of

1) mediate the earliest "rolling" of leukocytes on the lumenal surface of activated endothelium, and the interaction between leukocytes and activated platelets or other leukocytes found at sites of inflammation. CD162 is expressed on neutrophils, monocytes, and most lymphocytes including NK and T cells but PSGL-1 stains B cells at significantly lower levels than other cell types.

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com

Antigen 1. Snapp KR, *et al.* 1998. *Blood* 91:154. References:

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com