

Brilliant Violet 785™ anti-human HLA-A2

Catalog # / Size: 2316635 / 25 tests
2316640 / 100 tests

Clone: BB7.2

Isotype: Mouse IgG2b, κ

Reactivity: Human

Concentration: NULL

Applications:

Applications: Flow Cytometry

Recommended Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤ 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

Brilliant Violet 785™ excites at 405 nm and emits at 785 nm. The bandpass filter 780/60 nm is recommended for detection, although filter optimization may be required depending on other fluorophores used. Be sure to verify that your cytometer configuration and software setup are appropriate for detecting this channel. Refer to your instrument manual or manufacturer for support. Brilliant Violet 785™ is a trademark of Sirigen Group Ltd.

Application Notes: The BB7.2 antibody recognizes human leukocyte antigen (HLA) A2 which is a subset of MHC-class I molecules encoded by A*02 alleles.

Additional reported applications (for the relevant formats) include: immunoprecipitation³.

Description: HLA-A2 is most common in Northern Asia and North America populations. MHC class I antigens associated with $\beta 2$ -microglobulin are expressed by all human nucleated cells. MHC class I molecules are involved in presentation of antigens to CD8⁺ T cells, playing an important role in cell-mediated immune responses and tumor surveillance.