## PE anti-human CD89

	2370515 / 25 tests 2370520 / 100 tests	
Clone:	A59	A
Isotype:	Mouse IgG1, к	
Immunogen:	Ag8.653 myeloma cells	Number of the second se
<b>Reactivity:</b>	Human	ed ather cell Number
Preparation:	The antibody was purified by affinity chromatography and conjugated with PE under optimal conditions. The solution is free of unconjugated PE and unconjugated antibody.	$ \begin{array}{c} \hline aligned control \\ \hline aligned control $
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin USA).	Human peripheral blood granulocytes were stained with
Workshop Number:	V MR30	anti-human CD89 PE (clone A59) (filled histogram) or mouse IgG1, κ PE isotype control (open
<b>Concentration:</b>	Lot-specific	histogram).

## **Applications:**

 Applications: Flow Cytometry
 Recommended Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.
 Application References: 1. Monteiro RC, et al. 1992. J. Immunol. 148:1764.
 Shen L. 1992. J. Leukoc. Biol. 51:373.

 Schlossman S, et al. Eds. 1995. Leucocyte Typing V. Oxford University Press. New York.
 Rogers KA, et al. 2004. Immunology 113:178.

Antigen	1. Patry C, et al. 1996. J. Immunol. 156:4442.	
<b>References:</b>	2. de Wit, et al. 1995. J. Immunol. 155:1203.	
	3. Honorio-França AC, et al. 2001. J. Leukoc. Biol. 69:289.	

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com