## **Product Data Sheet**

## Brilliant Violet 510<sup>™</sup> anti-human CD56 (NCAM)

| Catalog # / Size:  | 2412665 / 25 tests<br>2412670 / 100 tests   | 4   |
|--------------------|---|---|
| Clone:             | 5.1H11  |   |
| Isotype:           | Mouse IgG1, κ   |   |
| Immunogen:         | Human myotube cells.  | z   |
| <b>Reactivity:</b> | Human   | he Cell   |
| Preparation:       | The antibody was purified by affinity<br>chromatography and conjugated with<br>Brilliant Violet 510 <sup>™</sup> under optimal<br>conditions. The solution is free of<br>unconjugated Brilliant Violet 510 <sup>™</sup> and<br>unconjugated antibody. | 0 10 <sup>2</sup> 10 <sup>3</sup> 10 <sup>4</sup> 10 <sup>5</sup><br>Log Fluorescence Intensity         |
| Formulation:       | Phosphate-buffered solution, pH 7.2,<br>containing 0.09% sodium azide and BSA<br>(origin USA).  | Human peripheral blood<br>lymphocytes were stained with<br>CD56 (clone 5.1H11) Brilliant                |
| Concentration:     | Lot-specific  | Violet <sup>™</sup> 510 (filled histogram) or<br>mouse IgG1, κ APC isotype control<br>(open histogram). |

## **Applications:**

| Applications:              | Flow Cytometry  |
|----------------------------|---|
| Recommended<br>Usage:      | Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is $\leq 5$ microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.   |
|                            | Brilliant Violet 510 <sup>™</sup> excites at 405 nm and emits at 510 nm. The bandpass filter 510/50 nm is recommended for detection, although filter optimization may be required depending on other fluorophores used. <b>Be sure to verify that your cytometer configuration and software setup are appropriate for detecting this channel.</b> Refer to your instrument manual or manufacturer for support. Brilliant Violet 510 <sup>™</sup> is a trademark of Sirigen Group Ltd.   |
| Application<br>Notes:      | For <i>in vivo</i> studies or highly sensitive assays, we recommend Ultra-LEAF <sup><math>m</math></sup> purified antibody (Cat. No. 362548) with a lower endotoxin limit than standard LEAF <sup><math>m</math></sup> purified antibodies (Endotoxin <0.01 EU/microg).   |
| Application<br>References: | 1. Walsh FS, <i>et al.</i> 1981. <i>Nature</i> 289:60. (FC)<br>2. Pavlath GK, <i>et al.</i> 1986. <i>J. Cell Biol.</i> 102:124. (FC)<br>3. Pavlath GK, <i>et al.</i> 1989. <i>Nature</i> 337:570. (FC)<br>4. Pulido R, <i>et al.</i> 1988. <i>J. Immunol.</i> 140:3851. (FC)  |
| Description:               | CD56 is a single transmembrane glycoprotein also known as NCAM (neural cell<br>adhesion molecule), Leu-19, or NKH1. It is a member of the lg superfamily. The<br>140 kD isoform is expressed on NK and NKT cells. CD56 is also expressed in the<br>brain (cerebellum and cortex) and at neuromuscular junctions. Certain large<br>granular lymphocyte (LGL) leukemias, small-cell lung carcinomas, neuronal<br>derived tumors, myelomas, and myeloid leukemias also express CD56. CD56<br>plays a role in homophilic and heterophilic adhesion via binding to itself or<br>heparan sulfate. |
|                            |   |

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com

| Antigen            | 1. Lanier L, <i>et al.</i> 1991. <i>J. Immunol.</i> 146:4421 |
|--------------------|--|
| <b>References:</b> | 2. Hemperly J, et al. 1990. J. Mol. Neurosci. 2:71           |
|                    | 3. Cremer H, <i>et al.</i> 1994. <i>Nature</i> 367:455.      |

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com