## Spark Violet™ 423 anti-human CD27

Catalog # / 356442 / 100 tests

Size: 356441 / 25 tests

Clone: M-T271

Isotype: Mouse IgG1, ĸ

Human T cells from a T-ALL patient. Immunogen:

Reactivity: Human

The antibody was purified by affinity **Preparation:** 

chromatography and conjugated with Spark Violet<sup>™</sup> 423 under optimal

conditions.

Formulation: Phosphate-buffered solution, pH 7.2.

containing 0.09% sodium azide and

BSA (origin USA)

Workshop **Number:**  V 5T CD27.03

**Concentration:** Lot-specific Human peripheral blood lymphocytes were stained with anti-human CD3 APC and antihuman CD27 (clone M-T271) Spark Violet™ 423 (left) or mouse IgG1, κ Spark Violet™ 423 isotype

control(right).



Applications: Flow Cytometry

Recommended

Usage:

Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 μL per million cells in 100 μL staining volume or 5 μL per 100 μL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

\* Spark Violet™ 423 has a maximum excitation of 400 nm and a maximum

emission of 415 nm.

**Application** 

Notes:

Additional reported applications (for the relevant formats) include: immunohistochemical staining of formalin-fixed paraffin-embedded frozen

tissue sections<sup>1</sup>, immunofluorescent staining<sup>2</sup>, and ELISA<sup>3</sup>.

**Application** References:

1. Ma S, et al. 2011. J. Virol. 85:165. (IHC)

2. Manzo A, et al. 2008. Arthritis Rheum. 11:3377. (IF)

3. Kato K, et al. 2007. Exp. Hematol. 35:434. (ELISA)

**Description:** 

CD27 is a 50-55 kD type I membrane protein also known as S152 and T14. It is a lymphocyte-specific member of the TNF-receptor superfamily. CD27 is expressed on medullary thymocytes, virtually all mature T cells, some B cells, and NK cells. CD27 binds to CD70, and plays a role in costimulation of T cell activation and regulation of B cell differentiation and proliferation. The cytoplasmic domains of CD27 have also been shown to interact with TRAF2 and TRAF5 to elicit NF-kB and SAPK/JNK activation.

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com

## **Antigen** References:

- 1. Knapp W, et al. 1989. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford University Press.
- 2. Schlossman S, et al. 1995. Leucocyte Typing V: White Cell Differentiation Antigens. Oxford University Press.
- 3. Hintzen R, et al. 1994. Immunol. Today 15:307.
- 4. Agematsu K, et al. 1995. J. Immunol. 154:3627.