Product Data Sheet

PerCP anti-mouse CD5

Catalog # / Size: $1103075 / 25 \mu g$

1103080 / 100 μg

Clone: 53-7.3

Isotype: Rat IgG2a, κ

Immunogen: Mouse thymus or spleen

Reactivity: Mouse

Preparation: The antibody was purified by affinity


chromatography, and conjugated with PerCP under optimal conditions. The solution is free of unconjugated PerCP

and unconjugated antibody.

Formulation: Phosphate-buffered solution, pH 7.2,

containing 0.09% sodium azide.

Concentration: 0.2



C57BL/6 mouse splenocytes stained with 53-7.3 PerCP

Applications:

Applications: Flow Cytometry

Recommended Usage:

Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of

this reagent is ≤ 0.25 microg per 10^6 cells in 100 microL volume. It is recommended that the reagent be titrated for optimal performance for each

application.

* PerCP has a maximum absorption of 482 nm and a maximum emission of 675

nm.

Application Notes:

Additional reported applications (for the relevant formats) include:

immunoprecipitation1, and immunohistochemistry2 of acetone-fixed frozen tissue

sections, zinc-fixed paraffin-embedded sections and formalin-fixed paraffin-

embedded sections.

Application References:

Ledbetter JA, et al. 1979. Immunol. Rev. 47:63. (IP)
Ledbetter JA, et al. 1980. J. Exp. Med. 152:280. (FC, IHC)

3. Bourdeau A. et al. 2007. Blood doi:10.1182/blood-2006-08-044370.

4. Plzakova L, et al. 2014. Microb Pathog. 75:49. PubMed

Description:

CD5 is a 67 kD protein, also known as Lyt-1, Ly-1, T1, Tp67, or Ly-12. It is a member of the scavenger receptor cysteine-rich protein superfamily (SRCR) and primarily expressed on thymocytes, T cells, and B-1 cells. Although mature α/β T cells express high levels of CD5, very few γ/δ T cells express this antigen. The interaction of CD5 with CD72, gp35-37, TCR, or BCR is involved in T and B cell activation.

Antigen References:

1. Barclay A, et al. 1997. The Leukocyte Antigen FactsBook Academic Press.

2. Kipps TJ. 1988. Adv. Immunol. 47:117.

3. Antin JH, et al. 1985. J. Immunol. 136:505.

4. Tarakhovsky A, et al. 1