
Product Data Sheet

APC/Cy7 anti-mouse CD14

Catalog # / Size:	1216585 / 25 μg 1216590 / 100 μg
Clone:	Sa14-2
Isotype:	Rat IgG2a, к
Immunogen:	Mouse thymus or spleen
Reactivity:	Mouse
Preparation:	The antibody was purified by affinity chromatography and conjugated with APC/Cy7 under optimal conditions. The solution is free of unconjugated APC/Cy7 and unconjugated antibody.
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
Concentration:	0.2




CD107b (Mac-3) FITC

Thioglycolate-elicited Balb/c mouse peritoneal macrophages were stained with CD107b (Mac-3) FITC and CD14 (clone Sa14-2) APC/Cy7 (top) or rat IgG2a, κ APC/Cy7 isotype control (bottom).

Applications:

Applications:	Flow Cytometry
Recommended Usage:	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤ 0.5 microg per million cells in 100 microL volume. It is recommended that the reagent be titrated for optimal performance for each application.



Description: CD14 is a 53-55 kD glycosylphosphatidylinositol (GPI)-linked membrane glycoprotein also known as LPS receptor. CD14 is expressed on macrophage dendritic cells, Kupffer cells, hepatocytes, and granulocytes. As a high-affir receptor for LPS-LBP (LPS-binding protein) complex, CD14, in association w Toll-like Receptor 4 (TLR4) or 2 (TLR2), is involved in the clearance of gram negative pathogens.
--

Antigen	1. Stocks S, <i>et al.</i> 1990. <i>Biochem. J.</i> 268:275.
References:	2. Akashi S, <i>et al.</i> 2003. <i>J. Exp. Med.</i> 198:1035.
	3. Matsuura K, <i>et al.</i> 1994. <i>J. Exp. Med.</i> 179:1671. 4. Liu S, <i>et al.</i> 1998.

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com