Purified anti-mouse TIGIT (Vstm3)

Catalog # / Size: 1310510 / 500 μg

1310505 / 50 μg

Clone: 1G9

Isotype: Mouse IgG1, κ

Reactivity: Mouse


Preparation: The antibody was purified by affinity

chromatography.

Formulation: Phosphate-buffered solution, pH 7.2,

containing 0.09% sodium azide.

Concentration: 0.5



Con A-stimulated (3 days) C57BL/6 mouse splenocytes were stained with CD3 APC and purified TIGIT (clone IG9, top) or mouse IgG1, κ isotype control (bottom), followed by anti-mouse IgG1 PE.


Applications:

Applications: Flow Cytometry

Recommended

Usage:

Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤1.0 microg per million cells in 100 microL volume. It is recommended that the reagent be titrated for optimal performance for each application.



Application References:

1. Joller N, et al. 2010. J. Immunol. 186:1338.

Description:

T cell immunoreceptor with Ig and ITIM domains (TIGIT), also known as V-set and transmembrane domain-containing protein 3 (Vstm3), is a 26 kD, type I transmembrane protein and member of the CD28 family. TIGIT is expressed on activated T cells, follicular T helper, memory, and regulatory T cells as well as on NK cells. Its binding partners include CD155 (PVR) and CD112 (PVRL2). TIGIT is a negative regulator of NK and T cell activation. Engagement of TIGIT by dendritic cells results in their differentiation into a tolerogenic phenotype, with an increased secretion of IL-10 and a diminished production of IL-12. Mice deficient for TIGIT are more susceptible to autoimmune disease.

Antigen References:

1. Levin SD, et al. 2011. Eur. J. Immunol. 41:902.

2. Yu X, et al. 2009. Nat. Immunol. 10:48.

3. Stanietsky N, et al. 2009. P. Natl. Acad. Sci. USA 106:17858.