

PerCP/Cyanine5.5 anti-mouse CD200R3

Catalog # / Size: 1311070 / 100 µg
1311065 / 25 µg

Clone: Ba13

Isotype: Rat IgG2a, κ


Immunogen: Mouse primary basophils

Reactivity: Mouse

Preparation: The antibody was purified by affinity chromatography and conjugated with PerCP/Cyanine5.5 under optimal conditions. The solution is free of unconjugated PerCP/Cyanine5.5 and unconjugated antibody.

Formulation: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.

Concentration: 0.2 mg/ml



BALB/c mouse peripheral blood leukocytes were stained with FcεRIα PE and CD200R3 (clone Ba13) PerCP/Cyanine5.5 (left) or Rat IgG2a, κ PerCP/Cyanine5.5 isotype control (right)

Applications:

Applications: Flow Cytometry

Recommended Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤ 0.25 µg per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application.

* PerCP/Cyanine5.5 has a maximum absorption of 482 nm and a maximum emission of 690 nm.

Application Notes: Ba13 recognizes circulating and bone marrow basophils; it also recognizes a subset of mast cells in the peritoneal cavity and skin. Additional reported applications (for the relevant formats) include: stimulation of bone marrow derived basophils to produce IL-4.

- Application References:**
1. Takakura N, et al. 1996. *J. Invest. Dermatol.* 107:770.
 2. Liao C, et al. 2010. *J. Clin. Invest.* 120:242. (Block)
 3. Chen H, et al. 2015. *ASN Neuro* 8:7. [PubMed](#)

Description: CD200R3, also known as CD200RLb and OX-2 Receptor 3, is a disulfide-linked dimeric CD200R-like receptor which belongs to immunoglobulin superfamily. Its positively charged amino acid lysine associates with ITAM- or YxxM motif-bearing adaptor molecules such as DAP12, DAP10, FcRγ, and CD3ζ. CD200R3 functions as an activating receptor to regulate IgE independent immune response.

- Antigen References:**
1. Voehringer D, et al. 2004. *J. Biol. Chem.* 52:54117.
 2. Kojima T, et al. 2007. *J. Immunol.* 179:7093.
 3. Sato K, et al. 2009. *Blood* 113:4780.