Product Data Sheet

PerCP/Cy5.5 anti-human CD45RA

Catalog # / Size:	2120605 / 25 tests 2120610 / 100 tests	
Clone:	HI100	al la
Isotype:	Mouse lgG2b, κ	ě M. 🕌
Reactivity:	Human	
Preparation:	The antibody was purified by affinity chromatography, and conjugated with PerCP/Cy5.5 under optimal conditions. The solution is free of unconjugated PerCP/Cy5.5 and unconjugated antibody.	10 ⁰ 10 ¹ 10 ² 10 ³ 10 ⁴
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin USA).	10 ⁰ 10 ¹ 10 ² 10 ³ 10 ⁴ Log Fluoresence Intensity Human peripheral blood lymphocytes stained with HI100
Workshop Number:	IV N906	PerCP/Cy5.5
Concentration:	Lot-specific	

Applications:

Applications:	Flow Cytometry
---------------	----------------

Recommended Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

* PerCP/Cy5.5 has a maximum absorption of 482 nm and a maximum emission of 690 nm.

Application Notes: Additional reported applications (for relevant formats of this clone) include: inhibition of CD45 functions2, immunohistochemical staining of frozen tissue sections3 and formalin-fixed paraffin-embedded tissue sections4, and immunofluorescence^{15,16}.

Application	1. Knapp W, <i>et al.</i> 1989. Leucocyte Typing IV. Oxford University Press. New York.
References:	2. Yamada T, <i>et al.</i> 2002. <i>J. Biol. Chem.</i> 277:28830. (WB, Block)
	3. Weninger W, <i>et al.</i> 2003 <i>J. Immunol.</i> 170:4638. (IHC)
	4. Imanguli MM, <i>et al.</i> 2009. <i>Blood</i> . 113:3620 (IHC)
	5. Roque S, <i>et al.</i> 2007 <i>. J. Immunol.</i> 178:8028. (FC) <u>PubMed</u>
	6. Smeltz RB. 2007. <i>J. Immunol.</i> 178:4786. (FC) PubMed
	7. Palendira U, <i>et al.</i> 2008. <i>Blood</i> (FC) <u>PubMed</u>
	8. Kuttruff S, <i>et al.</i> 2009. <i>Blood</i> 113:358. (FC) <u>PubMed</u>
	10. Thakral D, <i>et al.</i> 2008. <i>J. Immunol.</i> 180:7431. (FC) <u>PubMed</u>
	11. Alanio C, <i>et al.</i> 2010. <i>Blood</i> 115:3718. (FC) <u>PubMed</u>
	12. Iannello A, <i>et al.</i> 2010. <i>J. Immunol.</i> 184:114. (FC) <u>PubMed</u>
	13. Yoshino N, <i>et al.</i> 2000. <i>Exp. Anim. (Tokyo)</i> 49:97. (FC)
	14. Guereau-de-Arellan M, <i>et al.</i> 2011. <i>Brain</i> . 134:3578. <u>PubMed</u>
	15. Canque B, <i>et al.</i> 2000. <i>Blood</i> 96:3748. (IF)
	16. Imanguli MM, <i>et al.</i> 2009. <i>Blood</i> 13:3620. (IF)

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com **Description:** CD45RA is a 205-220 kD single chain type I glycoprotein. It is an exon 4 splice variant of the tyrosine phosphatase CD45. The CD45RA isoform is expressed on resting/naïve T cells, medullary thymocytes, B cells and monocytes. CD45RA enhances both T cell receptor and B cell receptor signaling. CD45 non-covalently associates with lymphocyte phosphatase-associated phosphoprotein (LPAP) on T and B lymphocytes. CD45 has been reported to be associated with several other cell surface antigens including CD1, CD2, CD3, and CD4. CD45 has also been reported to bind galectin-1. CD45 isoform expression can change in response to cytokines.

Antigen 1. Thomas M. 1989. Annu. Rev. Immunol. 7:339.

References: 2. Trowbridge I, et al. 1994. Annu. Rev. Immunol.12:85.

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com