## PerCP/Cy5.5 anti-human CD56 (NCAM)

Catalog # / Size: 2123130 / 100 tests

2123125 / 25 tests

Clone: MEM-188

Isotype: Mouse IgG2a, κ

KG-1 human acute myelogenous Immunogen:

leukemia cell line

Reactivity: Human

**Preparation:** The antibody was purified by affinity

chromatography and conjugated with PerCP/Cy5.5 under optimal conditions. The solution is free of unconjugated PerCP/Cy5.5 and unconjugated

antibody.


Formulation: Phosphate-buffered solution, pH 7.2,

containing 0.09% sodium azide and

0.2% (w/v) BSA (origin USA).

Workshop **Number:**  VI NK26

Concentration: Lot-specific



Human peripheral blood lymphocytes stained with MEM-188 PerCP/Cv5.5

## **Applications:**

**Applications:** Flow Cytometry

Recommended

Usage:

Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

\* PerCP/Cy5.5 has a maximum absorption of 482 nm and a maximum emission of

**Application** 

Notes:

Additional reported applications (for the relevant formats) include:

immunoprecipitation, immunohistochemical staining of formalin-fixed paraffin-

embedded tissue sections, and Western blotting (non-reducing).

**Application** 

1. Kishimoto T, et al. Eds. 1997. Leucocyte Typing VI. Garland Publishing Inc.

References: London.

**Description:** CD56 is a single transmembrane glycoprotein also known as N-CAM (Neural Cell

Adhesion Molecule), Leu-19, or NKH1. It is a member of the Ig superfamily. The 140 kD isoform is expressed on NK cells and NK-T cells. CD56 is also expressed in the brain (cerebellum and cortex) and at neuromuscular junctions. Certain large granular lymphocyte (LGL) leukemias, small-cell lung carcinomas, neuronal derived tumors, myelomas, and myeloid leukemias also express CD56. CD56 plays a role in homophilic and heterophilic adhesion via binding to itself or

heparin sulfate.

**Antigen** References: 1. Lanier L, et al. 1991. J. Immunol. 146:4421.

2. Hemperly J, et al. 1990. J. Mol. Neurosci. 2:71.

3. Cremer H, et al. 1994. Nature 367:455.

