## PerCP/Cy5.5 anti-human CD10

-	2161080 / 100 tests 2161075 / 25 tests	□ Human pre-B cell line, Nalm-6, stained with HI10a PerCP/Cy5.5
Clone:	HI10a	
lsotype:	Mouse IgG1, к	
<b>Reactivity:</b>	Human	
Preparation:	The antibody was purified by affinity chromatography and conjugated with PerCP/Cyanine5.5 under optimal conditions. The solution is free of unconjugated PerCP/Cyanine5.5 and unconjugated antibody.	
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin USA).	
Workshop Number:	V CD10.7	
<b>Concentration:</b>	Lot-specific	

## **Applications:**

Applications:	Flow Cytometry	
Recommended Usage:	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 $\mu$ l per million cells or 5 $\mu$ l per 100 $\mu$ l of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.	
	* PerCP/Cyanine5.5 has a maximum absorption of 482 nm and a maximum emission of 690 nm.	
Application Notes:	Additional reported (for the relevant formats) applications include: immunohistochemistry <sup>6</sup> .	
Application References:	1. Shipp M, <i>et al.</i> 1993. <i>Blood</i> 82:1052. 2. Lu B, <i>et al.</i> 1995. <i>J. Exp. Med.</i> 181:2271.	
Description:	CD10 is a 100 kD neutral endopeptidase and a member of the metalloprotease family. It is a type II transmembrane protein also known as common acute lymphoblastic leukemia antigen (CALLA), enkephalinase, and neprilysin. CD10 is expressed on B cell precursors, T cell precursors, and neutrophils. CD10 is involved in B cell development and has been shown to	

biologically active peptides.Antigen1. Shipp M, et al. 1993. Blood 82:1052.References:2. Lu B, et al. 1995. J. Exp. Med. 181:2271.

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com

bind opioid enkephalins, bradykinin, angiotensins I and II, and other