Product Data Sheet

PE/Dazzle™ 594 anti-human CD21

Catalog # / 2374605 / 25 tests

Size: 2374610 / 100 tests

Clone: Bu32

Isotype: Mouse IgG1, κ

Reactivity: Human, Non-human primate

Preparation: The antibody was purified by affinity

chromatography and conjugated with PE/Dazzleâ, \$ 594 under optimal conditions. The solution is free of unconjugated PE/Dazzleâ, \$ 594 and

unconjugated antibody.


Formulation: Phosphate-buffered solution, pH 7.2,

containing 0.09% sodium azide and

0.2% (w/v) BSA (origin USA).

Workshop Number: V CD21.4, VI CD21.5

Concentration: Lot-specific



Human peripheral blood lymphocytes were stained with CD19 BV421 and CD21 (clone Bu32, left) PE/Dazzle™ 594 or Mouse IgG1, κ PE/Dazzle™ 594 isotype control (right).

Applications:

Applications: Flow Cytometry

Recommended

Usage:

Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 μ l per million cells or 5 μ l per 100 μ l of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

* PE/Dazzle™ 594 has a maximum excitation of 566 nm and a maximum

emission of 610 nm.

Application Notes:

Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen sections 4 .

Application References:

1. Kishimoto T, Eds. 1997. Leukocyte Typing VI. Garland Publishing Inc.

2. Moore MD, et al. 1987. Proc. Natl. Acad. Sci. USA 84:9194.

3. Szakonyi G, et al. 2001. Science 292:1725.

4. Weis JJ, et

Description: CD21 is a 145 kD transmembrane protein also known as complement C3d

receptor (C3dR), complement receptor 2 (CR2), and Epstein-Barr virus receptor. CD21 is expressed on B cells, follicular dendritic cells, subsets of normal thymocytes and T cells, and some epithelial cells. CD21 is the receptor used by Epstein-Barr virus to infect B cells and is also the

complement receptor for C3d. CD21 has also been shown to interact with a number of proteins, including CD23, CD19, annexin VI, CD81, iC3b,

complement receptor 1 (CR1, CD35), and interferon-alpha 1 (IFN- α 1).

Antigen References:

1. Kishimoto T, Eds. 1997. Leukocyte Typing VI. Garland Publishing Inc.

2. Moore MD, et al. 1987. Proc. Natl. Acad. Sci. USA 84:9194.

3. Szakonyi G, et al. 2001. Science 292:1725.

4. Weis JJ, et al. 1984. Proc. Natl. Acad. Sci. USA 81:881.