Product Data Sheet

PerCP/Cy5.5 anti-human CD33

Catalog # / Size:	2433080 / 100 tests 2433075 / 25 tests	100-
Clone:	P67.6	
Isotype:	Mouse lgG1, к	
Immunogen:	FMY9S5 cells expressing CD33 gene.	000 000 000 000 000 000 000 000 000 00
Reactivity:	Human	aditive Cell Number
Preparation:	The antibody was purified by affinity chromatography and conjugated with PerCP/Cy5.5 under optimal conditions. The solution is free of unconjugated PerCP/Cy5.5 and unconjugated antibody.	20 0 10 ² 10 ³ 10 ⁴ 10 ⁵ Log Fluorescence Intensity
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin USA).	Human peripheral blood monocytes were stained with CD33 (clone P67.6) PerCP/Cy5.5 (filled histogram) or mouse IgG1, κ PerCP/Cy5.5 isotype control (open histogram).
Concentration:	Lot-specific	

Applications:

Recommended Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 microL per million cells or 5 microL per 100 microL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

* PerCP/Cy5.5 has a maximum absorption of 482 nm and a maximum emission of 690 nm.

Application 1. Hoyer J, *et al.* 2008. *Am. J. Clin. Pathol.* 129:316. References:

Description: CD33, also known as Siglec-3, gp67, and p67, is a 67 kD type I transmembrane glycoprotein. It is a sialoadhesion immunoglobulin superfamily member, which is expressed on myeloid progenitors, monocytes, granulocytes, dendritic cells, and mast cells. CD33 is absent on normal platelets, lymphocytes, erythrocytes, and hematopoietic stem cells. CD33 functions as a sialic acid-dependent cell adhesion molecule with carbohydrate/lectin binding activity.

Antigen	1. Favaloro E, <i>et al.</i> 1988. <i>Br. J. Haematol.</i> 69:163.
References:	2. Freeman S, <i>et al.</i> 1995. <i>Blood</i> 85:2005.

For research use only. Not for diagnostic use. Not for resale. Sony Biotechnology Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products. Sony Biotechnology Inc. 1730 North First Street, San Jose, CA 95112 www.sonybiotechnology.com